

SC G&M 2000 SRL
Bucuresti

Utilizarea dezinfectantului CLOROM

Dezinfectantul CLOROM se prezinta sub forma unor tablete. Este dezinfectant de uz general. Solutiile de lucru se obtin prin dizolvarea tabletelor in apa.

Concentratie	Actiune	Timp contact, minute
1 tb./3 L apa	Bactericid (inclusiv MRSA) fungicid si levuricid	5
1 tb./5 L apa	Bactericid si levuricid	5
1 tb./3 L apa	Bactericid, levuricid, Fungicid, virucid, sporicid (inclusiv Clostridium difficile), tuberculocid si micobactericid	60
1 tb./3 L apa	Virucid	20

ADMINISTRATOR

ISTRATE Gheorghe



INGINER CHIMIST

Pascu Ecaterina

